

TECHNICAL DATA TREATMENT LIGHTS


FITTER WITH	4 1501 - 000070000
FITTED WITH	1 x LED lamp - D68873000
LUMINOUS POWER	55 000 Lux / 1,0 m
	55 000 Lux / 1,0 m
ILLUMINANCEFIELD 1,0 M	Ø 14-32 cm
COLOR RENDERING INDEX (CRI)	> 93
COLOUR TEMPERATUR	approx. 3700 K / 4500 K / 5000 K
DEPTH OF ILLUMINATION	1,074 m
L1+L2	
WORK EQUIPMENT	electronic converter
CONNECTED LOAD	100-240 V; 50/60 Hz
MAINS LEAD	direct power supply
LUMINAIRE BODY	aluminium
WEIGHT (NET)	approx. 25.0 kg
POWER CONSUMPTION	approx. 30 W
BALANCE OF ARTICULATION JOINTS	spring (weight compensation)
CLASS OF PROTECTION	I
ACCESSORY	sterilizable handle, adapter for suspended ceiling
NORMS	EN 60601-1, EN 60601-2-41
SPECIAL FEATURES	3 color shades, 6 dimming levels, sterilisable grip, variable illumination field size, maintenance free
FASTENING	ceiling mount, when mounting ceiling tube could be cut to length 30-120 cm, intermediate ceiling adapter (on request)
ORDER NUMBER	D15071000 - pure white





LIGHT DISTRIBUTION TECHNICAL DRAWING











illuminancefield 14 cm illuminancefield 32 cm

Photometric data -10% / +20% tolerance