

# Dmed® saturn LED 30 W

## TECHNICAL DATA

## TREATMENT LIGHTS

FITTED WITH	1 x LED lamp - D68873000
LUMINOUS POWER	55 000 Lux / 1,0 m 55 000 Lux / 1,0 m
ILLUMINANCEFIELD 1,0 M	Ø 14-32 cm
COLOR RENDERING INDEX (CRI)	> 93
COLOUR TEMPERATUR	approx. 3700 K / 4500 K / 5000 K
DEPTH OF ILLUMINATION L1+L2	1,074 m
WORK EQUIPMENT	electronic converter
CONNECTED LOAD	100-240 V; 50/60 Hz
MAINS LEAD	direct power supply
LUMINAIRE BODY	aluminium
WEIGHT (NET)	approx. 18.0 kg
POWER CONSUMPTION	approx. 30 W
BALANCE OF ARTICULATION JOINTS	spring (weight compensation)
CLASS OF PROTECTION	I
NORMS	EN 60601-1, EN 60601-2-41
SPECIAL FEATURES	3 color shades, 6 dimming levels, sterilisable grip, variable illumination field size, maintenance free
FASTENING	wall mount
ORDER NUMBER	D15133000 - pure white



## LIGHT DISTRIBUTION

## TECHNICAL DRAWING

illuminancefield 14 cm

illuminancefield 32 cm

