

# Dmed® visano 20-2 C T1

## TECHNICAL DATA

## EXAMINATION LIGHTS

FITTED WITH	1 x LED lamp
LUMINOUS POWER	60 000 Lux / 0.5 m
ILLUMINANCEFIELD M	Ø 19 cm
COLOR RENDERING INDEX (CRI)	> 93
COLOUR TEMPERATUR	approx. 3500 K, 4500 K
WORK EQUIPMENT	electronic converter
CONNECTED LOAD	90-277V; 50/60Hz
MAINS LEAD	direct power supply
WEIGHT (NET)	approx. 7.5 kg
POWER CONSUMPTION	approx. 20 W
USAGE	touch panel
CLASS OF PROTECTION	II
ACCESSORY	adapter for suspended ceiling
NORMS	EN 60601-1, EN 60601-2-41
SPECIAL FEATURES	2 color shades, 4 dimming levels, antimicrobial, Rotary feedthrough (360° Turnover axis)
FASTENING	ceiling mount
SYSTEM OF PROTECTION	IP 43
ORDER NUMBER	D15450000 - pure white



## LIGHT DISTRIBUTION

## TECHNICAL DRAWING

