

BELLA *energy efficient and ecologic*

Lamp perfectly designed for minor surgery, dermatology and intensive cares...



- Lamp's power 17 W
- 5500° Kelvins
- 1440 Lumens
- Guarantee : 3 years
- LEDs' lifetime 50 000 hours (No UV output and not much infrared)
- Color Rendering Index (IRC) 95
- Cable length 3 meters allowing maximum mobility
- Insulation Class II
- Available colors: white
- Lamp's head dimensions 160 x 133 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Distance 30 cm: 55000 lux, Illuminated field: 8 cm
- Distance 50 cm: 42000 lux, Illuminated field: 10 cm
- Distance 100 cm: 13000 lux, Illuminated field: 20 cm
- Flexible length: 700 mm
- Double intensity switch in option (on the top of the lamp)
- Switch on/off or no contact switch
- The ergonomic handle allows a precise and intuitive positioning.
- Easy cleaning and disinfection thanks to the smooth surfaces and quality materials
- LED are implanted on a printed circuit cooled by an aluminum heat sink and not by a fan which is cause of unreliability, noise, and attracts dust.
- Codes : LED17650 – LED17650DI



lid-lamps.com