

CARLA *energy efficient and ecologic*



Examination lamp for localized and specific examination such as gynecology and dermatology
(Optic angle 25°)

- 4,2 watts
- 5500° kelvins
- 360 Lumens
- Guarantee 3 years
- LEDs' lifetime 50 000 hours (Cold white light, no UV or infrared output)
- Color Rendering Index (IRC) 95
- Cable length 3 meters allowing maximum mobility
- Insulation Class II
- Available colors: white, black
- Lamp's head dimensions 60 x 95 mm
- Lamp's base dimensions 69 x 69 mm, H: 60 mm, inter axial: 54,5 mm
- Distance 30 cm: 8000 lux, illuminated field: 20 cm
- Distance 50 cm: 3300 lux, illuminated field: 25 cm
- Distance 100 cm: 800 lux, illuminated field: 50 cm
- Flexible length: 650 mm or 1000 mm
- Switch on/off or no contact switch
- Option with battery: 2h autonomy (charger supplied)
- Code : LEDB0650 – LEDB0900 – LEDB0650-ISC
LEDB0900-ISC – LEDBAT650 – LEDBAT900



lid-lamps.com