

HEPTA *energy efficient and ecologic*

Very polyvalent examination lamp, powerful and comfortable, without any heat output, optimum handling



- 7 watts
- 5500° kelvins
- 840 Lumens
- Guarantee 3 years
- LEDs' lifetime 50 000 hours (Cold white light, no UV or infrared output)
- Color Rendering Index (IRC) 95
- Cable length 3 meters allowing maximum mobility
- Insulation Class II
- Available colors: white, black
- Lamp's head dimensions 110 x 110 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Distance 30 cm: 47500 lux, illuminated field: 8 cm
- Distance 50 cm: 25000 lux, illuminated field: 10 cm
- Distance 100 cm: 7000 lux, illuminated field: 20 cm
- Flexible length: 650 mm or 1000 mm
- Switch on/off or dimmer or no contact switch
- Option with battery: 2h autonomy (charger supplied)
- Code : LED07651-LED07651V-LED071001-LED071001V
LED06BAT650-LED06BAT1000



lid.lamps.com